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(EN) Instruction Manual

Dewaxing Chamber Furnaces

CONTENTS		Page
1.	PREFACE	3
2.	IMPORTANT SAFETY INSTRUCTIONS	3
2.1.	General information	3
2.2.	General safety instructions	3
2.3.	Technical data	4
2.4.	Protective clothing and safety equipment	5
2.5.	Operating safety instructions	5
3.	START-UP	6
3.1.	Delivery	6
3.2.	Installation environment / Location	6
3.3.	Assembly of frontloader dewaxing kiln	6
3.4.	Connecting the frontloader dewaxing kiln	7
3.5.	Power supply	7
3.6.	Frontloader dewaxing kiln initial firing	7
4.	GENERAL OPERATING INSTRUCTIONS	8
4.1.	Operating elements	8
4.2.	Description of kiln unit	8
4.3.	Loading the kiln	10
4.4.	Firing - starting to operate the kiln unit	11
4.5.	Correct operation during firing	11
5.	THE WAX-FUSION PROCESS	12
6.	BURNING OUT PROCESS	15
7.	MAINTENANCE / CARE and CLEANING	16
7.1.	Maintenance frontloader dewaxing kiln	16
7.2.	Replacing the heating elements	16
7.3.	Care and cleaning	17
8.	TROUBLESHOOTING TIPS	17
8.1.	Firing chamber	17
8.2.	Kiln control	17
9.	WARRANTY OBLIGATIONS	18
10.	PROPERTY RIGHTS / TRADE NAMES / DISCLAIMER	19
11.	CONTACTS / ASSISTANCE	19

1. PREFACE

Congratulations, you have chosen a dewaxing chamber kiln - a high-quality product, meeting highest requirements.

This instruction manual will help you to familiarise yourself with your new kiln. We have put together some important information and guidelines that will make operating your kiln as safe and simple as possible.

Please read the instruction manual carefully before using your kiln.

Make sure you understand the features and functions of your kiln and control unit.

2. IMPORTANT SAFETY INSTRUCTIONS

2.1. General information

Please make sure that you fully understand both the safety instructions and the safety icons, in order to eliminate potential dangers. Read the following safety instructions carefully before starting to operate your frontloader dewaxing kiln.

Keep your instruction manual available at all times. For your own safety only use original spare parts!

2.2. General safety instructions

CAUTION:
Hot surface!

Caution: Electrically live
components - Danger to Life!

Caution: Disconnect
power plug before
opening the switch box!

2.3. Technical data

Dewaxing Chamber Furnaces **KE-EW** series to 1000°C

Model	T _{max}	Internal dimensions mm			External dimensions mm			Output	Supply	Connector	Voltage	Weight
Volume	°C	b	t	h	B	T	H	kW	A	Plug	V	kg
KE 6/10 EW	1000	180	170	220	380	550	780	1.5	6.5	Schuko	230 V/50 Hz	45
KE 14/10 EW	1000	250	250	250	450	600	800	2.3	10.0	Schuko	230 V/50 Hz	55
KE 20/10 EW	1000	300	300	200	530	670	750	3.0	13.0	Schuko	230 V/50 Hz	60
KE 50/10 EW	1000	300	450	315	520	840	865	5.0	11.0	CEE 16 A	400 V/50 Hz	90
KE 150/10 EW	1000	460	470	580	800	1100	1800	13.5	20.0	CEE 32 A	400 V/50 Hz	460
KE 330/10 EW	1000	590	720	700	1000	1350	2100	22.0	32.0	CEE 32 A	400 V/50 Hz	680
KE 600/10 EW	1000	710	870	920	1100	1500	2300	40.0	59.0	CEE 63 A	400 V/50 Hz	1000

- Unit:** Electric frontloader dewaxing kiln EW series
- Moulded collecting tank on firing chamber floor
- Moulded collecting drawer below the kiln
- Controller:** Controller TC2088e-S-1
- Program memory for 15 programs
- Each program with 15 segments
- Clearly structured and easy-to-read display
- Heating:** 2-side heating (KE 6/10 EW – KE 50/10 EW)
4-side heating (from KE 150/10 EW)
- Maximum temperature:** 1000° C
Operating temperature: 950° C
- Warranty:** 24-months warranty (not on parts subject to wear)

Environmental conditions:

- Maximum environmental temperature: -5°C to +45°C
- Relative humidity: <80%, non-condensing
- Firing cycle: cold/cold
- Firing chamber atmosphere: oxidising

Unit designed and manufactured in accordance with DIN EN 746-1

Intended Use:

The KE-EW kiln can be used to burn out wax and plastic molded parts. All other use, e.g. the heating of food, the burning of other materials, etc., is considered as incorrect. The manufacturer accepts no liability for personal injury or property damage.

2.4. Protective clothing and safety equipment

Operating personnel must wear suitable protective clothing when opening the dewaxing kiln.

High temperatures are released that might cause physical injury and material damage.

Helmut Rohde GmbH does not accept liability for personal injury and material damage!

2.5. Operating safety instructions

The ROHDE frontloader dewaxing kiln can only be operated safely if you follow these safety instructions:

CAUTION:

When operated in an industrial unit, frontloader dewaxing kiln and controller must undergo a safety check to ensure correct functionality. This should be carried out by a qualified electrician before the initial operation and then at intervals in accordance with „DGUV Vorschrift 3“.

CAUTION:

The system and the burn-out process must be integrated by the customer into the safety concept of the operator. An analysis of the hazards at the workplace and in the work process has to be prepared. The safety instructions must be observed!

CAUTION:

Maintenance and repair of electronic components must be carried out by a qualified electrician.

CAUTION:

For safety reasons the kiln must be disconnected from the mains supply before any maintenance work is carried out.

CAUTION:

Never exceed the maximum temperature of 1000° C.

CAUTION:

The kiln must not be operated with an extension cable!

CAUTION:

When firing materials which release hazardous gases and fumes, an exhaust air system must be installed that directs these into the open air.

CAUTION:

The operator of the oven is instructed to comply with the applicable emission protection regulations (TA-Luft) "First General Administrative Regulation to the Federal Immission Control Act" of the German Federal Government! Or the applicable national laws!

3. START-UP

3.1. Delivery

The ROHDE frontloader dewaxing kiln will usually be delivered by a freight-forwarding agent. The goods must be checked immediately after delivery for visible damage on the packaging. Should you detect any damage, unpack the pallet together with the driver and recheck the goods again for damage. If you detect any damage please enter details on the delivery note and let the driver countersign your remarks. Keep one copy of the complaint for yourself. Inform the freight-forwarding agency immediately of the damage. Complaints submitted at a later date cannot be taken into consideration.

3.2. Installation environment / Location

When selecting a suitable place for your kiln, please note the following guidelines and prepare the kiln environment accordingly:


General instructions:

- Place the frontloader dewaxing kiln on an even surface.
- The distance to the walls should be at least 50 cm on each side.
- The floor, ceiling insulation, walls, dividing walls, panelling, etc. must be made of flame resistant material.
- Make sure that the kiln environment can be properly ventilated. If this is not the case, a ventilation system must be installed. Please consult a qualified ventilation specialist to find out whether a ventilation system is necessary.

3.3. Frontloader dewaxing kiln assembly

The frontloader dewaxing kiln is connected and put into operation by service technicians provided by the customer.

Make sure that the kiln is placed in a stable position on the workshop table. If required, the kiln legs can be readjusted.

Please loosen the screw connection. Now select the correct height setting for a stable kiln position and tighten the screw connection.

3.4. Connecting the frontloader dewaxing kiln

The frontloader dewaxing kiln is connected by the customer.

3.5. Power supply

The frontloader dewaxing kiln is connected by service technicians to the power supply provided by the customer.

Instructions Power connection / Residual current protective device (RCD)

If the kiln is to be operated in workshops or laboratories, a separate power supply with fuse protection must be installed by a qualified electrician.

Residual current protective devices (RCD) carrying a tripping current of 0.03 A (such as that used in damp rooms in flats) tend to trip early due to the high humidity of the rooms or fired goods.

A larger sized RCD can be selected (we recommend 0.3 A) provided that the respective circuit is used only for the kiln. If this cannot be guaranteed, a fixed power connection must be provided.

3.6. Frontloader dewaxing kiln initial firing

Carry out a dry firing before starting to use the kiln in in every day production. The "burning-in" by means of a dry firing is important in order to eliminate residual moisture from the kiln walls. It also generates a protective oxide layer on the heating elements which will considerably enhance the service life of these components.

Firing curve for initial firing:

- Heat up at 100° C/h
- End temperature 1000°C
- Holding time 1 h 30 min.

4. GENERAL OPERATING INSTRUCTIONS


4.1. Operating elements

Controller TC2088e-S-1

Please read the instruction manual TC 2088e in the second part of this documentation.


4.2. Description of kiln unit


Discharge grating inside the kiln


Clean the discharge grating after each heat treatment process.


Moulded collecting drawer inside the kiln


Clean the collecting drawer after each heat treatment process.


Kiln door hinge. Carefully pull the handle up and open the door. When closing the kiln door pull the handle down until it locks..

Exhaust air opening with exhaust air handle and exhaust air socket. Pull the handle to the left to close the exhaust air opening. Pull the handle to the right to open the exhaust air opening.

Air supply handle, lockable

Screw the exhaust air socket (marked red) into the hole on the ceiling of the kiln. The opening (see arrow) is located in a position that will allow fumes and gases to be released through an exhaust air socket (optional accessory). Plug the exhaust air tube into the exhaust air socket and use the fixing screw to fix it to the socket.

4.3. Loading the kiln

CAUTION:

The kiln may only be loaded by trained personnel! Operating personnel must wear suitable protective clothing and shoes. Use suitable kiln furniture when loading the kiln.

The manufacturer does not accept liability for personal injury and material damage!

CAUTION:

Hot surface! Danger of burning!

CAUTION:

The melting kiln and all its electrical components may only be operated by appropriately trained personnel!

The operating personnel must wear suitable protective clothing!
The manufacturer does not accept liability for personal injury and material damage!

CAUTION:

For the hot loading process (also hot unload), the main power switch must be switched OFF.

After closing the door and restarting the controller, the controller heats up again to the holding temperature (discharge temperature).


4.4. Firing - starting to operate the kiln unit

Step 1

Connect the kiln to the mains supply

Step 2

Fully close the kiln door

Step 3

Turn the switch key to position "I"
= ON

You can now start the firing from the kiln control TC 2088

4.5. Correct operation during firing

CAUTION:
Hot surface! Danger of burning!

- Do not place flammable objects near the kiln.
- The melting kiln may only be used in well ventilated areas.
- The melting kiln must be placed in a free-standing position in the room. Make sure that the heat release is not blocked. Do not place any objects on top of, or around, the kiln.
- When firing materials which release hazardous gases and fumes, an exhaust air system must be installed that directs these into the open air.
- Never use your kiln for firing inflammable materials or food.
- Only the operator may be present within the operating area. Unauthorised persons are not allowed.

5. THE WAX-FUSION PROCESS

Recommendations for the wax fusion process with the Rohde KE -EW

Basically, the process is divided into 7 steps:

Step: 1	Drying the cuvette
Step: 2	Melting out of the molding wax
Step: 3	Burning out the wax remains
Step: 4	Ready-burning or hardening of the mold
Step: 5	Hold on to finished firing temperature
Step: 6	Cooling to pouring or removal temperature
Step: 7	Keep at casting temperature

The next process step is ideally carried out after (actual) completion of the previous step. The thermal, physical and chemical processes take a certain amount of time and energy.

This means that the decisive factor for success is the process in the cuvette!

For the individual process steps, the following processes and indexing conditions result:

Step 1:

Drying of free water in the cuvette, a significant increase of the temperature will only occur in the dry parts of the cuvette. Too fast drying or drying at high temperatures (well above the boiling point of the water) should be avoided as the mold will otherwise be "overcooked". The cuvette is dry if, while drying and under consideration that no wax is running out, no loss of weight is detectable.

It is also possible to measure the temperature of the cuvette at a place with a large wall thickness, when it is equal to the furnace chamber temperature, the drying process is completed.

Step 1 & 2 merge into each other, whereby only a completely dried cuvette ensures a complete melting out of the wax.

Step 2:

The **melting out** of the wax should take place as continuously as possible from the sprue. The smelting process is completed when no more wax runs out. Here a visual inspection is possible.

Step 3:

The **burning** out serves to remove the wax residues from the pores of the molding material into which it has penetrated during melting. Especially forms with undercuts still significant amounts of wax must be burned. The oxygen required for this purpose must be supplied to the furnace by means of purging air, the reaction products being discharged via the exhaust gas flap / exhaust air flap.

The mold waxes used are mixtures of substances whose constituents have different evaporation and ignition temperatures, this being beneficial for a uniform course of the process step if not required time and sufficient scavenging air. If heating is too fast due to lack of air, cooking and soot formation in the furnace as well as cracking in the cuvettes due to thermally induced stresses occur.

In our experience, a flat ramp is more advantageous than holding at a constant temperature. An actual burnout is recognizable by an increase in the temperature of the cuvette above the controller setpoint, too rapid a burnout due to a sudden rise in the oven cavity temperature.

The end of the process step is reached when the cuvette temperature falls below the regulator set point again or no more combustion products over the chimney.

Step 4:

For **final firing** or hardening of the mold, the scavenging air quantity can be reduced, the exhaust air flap can initially remain open for the removal of the exhaust gases of the residual combustion. In the higher temperature range it should be closed for energy saving. The heating ramp should be chosen so that the stocking can also follow, to protect the heating elements of the furnace full load curves should be avoided.

Step 5:

The **holding** serves to equalize the temperature in the furnace chamber and the cuvettes and gives the running ceramic and mineralogical processes sufficient energy input and time.

Step 6:

The **cooling** to the removal temperature should not be too fast to avoid cracking in the molds. In a forced cooling with fan support is to be considered that the temperature of the mold core can still be significantly higher than the furnace chamber temperature. Again, the measurement of the cuvette temperature in the trial phase makes sense.

Step 7:

For **removal** when holding at casting temperature, the exhaust air flap can be closed again to save energy. A renewed rise in temperature above the holding temperature is an indication of a previous, too steep cooling ramp.

Before removal, a certain amount of time should be allowed for temperature compensation.

The user now has the problem of converting these steps into a brewing recipe, which is best achieved by observing the processes during a trial phase.

Here, the temperatures and times for the individual process steps must be determined to ensure a safe move on to the next step.

Different molding compounds and mold waxes (which we do not know!) Show a different behavior, a recipe that works with a cuvette in a stocking with two cuvettes of the same type (e.g double stocking: double mass to be heated and double mass to be melted out, if necessary double burning amount of wax) lead to significant problems. In particular, a sagging of the cuvette temperature and thus a shift of the running processes into the higher temperature range with an excessively programmed firing curve is fatal at the latest when the wax burns out.

So start the testing with little stocking, shapes and molding compounds with which you already have experience, simple, well-behaved shapes and enough safety margins at heating rates and (holding) times. Only make one change to the process and document input parameters and results.

Watch the process!

Please observe whether the oven and the stocker are following the programmed curve and, more important, whether the process steps are proceeding as intended. Divide the process into individual steps and take samples!

If, for example, a cuvette is not (yet) optimally melted out, it makes no sense to test parameter for the burnout, since the wax quantity, air requirement and possible heating rate will still change! On the other hand, well-melted cuvettes may enable you to increase the gross wax volume and thus the number of cuvettes or moldings in the kiln.

First, the process must work in the basics

An adaptation to the production requirements such as utilization of the maximum stocking weight / the maximum amount of wax, optimization of the time sequence and, if necessary, postponement of the critical process phases to the night hours can then take place systematically in the second step due to the knowledge gained.

6. BURNING OUT PROCESS

Polystyrene is melted out of prospective moulds made from plaster or similar material. The melted-out material is released into the discharge grating on the kiln floor and discharged into the collecting drawer by a drain. Polystyrene burns in this kiln atmosphere.

CAUTION:

During the process of Polystyrene melting inside the kiln, the temperature should not significantly exceed the melting point. There is a risk of igniting large amounts of Polystyrene which would cause damage to the kiln lining and the heating elements, This results from carbon which has not been completely burnt.

The manufacturer does not accept liability for personal injury and material damage!

Please contact your raw material supplier regarding the burning point.
The moulds can be burnt out after the Polystyrene has been burnt completely up to a maximum temperature of 900°C.

CAUTION:

It is essential to supply sufficient oxygen inside the kiln during the whole burning out process in order to guarantee complete combustion. There is a risk of deflagration!

The manufacturer does not accept liability for personal injury and material damage!

Use the designated air supply and exhaust air openings to ventilate the firing chamber.

The ventilation system to be installed must allow for an exhaust air volume of 30m³/h and may not at all impede the ventilation of the kiln.

The diagram below shows a burning out curve illustrating the above information.

7. MAINTENANCE / CARE and CLEANING

7.1. Maintenance frontloader dewaxing kiln

Maintenance operations:

CAUTION:

Electrically live components - Danger to Life! The kiln must be disconnected from all power supplies!

7.2. Replacing the heating elements

CAUTION:

Electrically live components - Danger to Life! The kiln must be disconnected from all power supplies!

Heating elements are subject to wear. Their resistance (Ohm) increases with each firing. Over the course of time this will lead to delays in the firing cycle due to a drop in performance, especially in the upper temperature range. If there is excessive wear we recommend that you replace the complete set of heating elements rather than just single elements. Replacing individual elements might lead to variations in temperature inside the kiln.

Caution: Heating elements which have already undergone firing are brittle and break easily (which is exempt from warranty). Therefore do not make any adjustments to the heating elements!

Useful information:

The service life of heating elements is influenced by:

- The kiln atmosphere: for example aggressive materials, particularly solvent-based materials or materials containing large amounts of sulphur or fluorine. It is important to ventilate the firing chamber properly.
- A tip for the firing professional: Always keep a spare set of heating elements! Should an emergency arise, this will save you unnecessary delay and allow you to continue firing as quickly as possible.

7.3. Care and cleaning

CAUTION:

We recommend an empty firing after every 5th firing. This will "clean" the heating elements allowing the oxide layer to renew itself. This will extend the service life of the elements

Please make sure that no impurities get into the firing chamber or come into contact with the heating elements. This will cause the heating elements to malfunction during subsequent firings. If, however, impurities do get into the heating chamber or onto the heating elements, clean them immediately, as burned-in impurities, etc. will damage the heating elements and bricks. If there is substantial damage, please contact Helmut Rohde GmbH or your retailer.

Clean the kiln regularly using a vacuum cleaner. This will also increase the service life of your heating elements.

Avoid reduction firing, as this will cause the oxidation layer to decompose, thus significantly reducing the service life of the heating elements.

8. TROUBLESHOOTING TIPS

8.1. Firing chamber: the kiln does not heat up at all or only very slowly

The kiln door must be fully closed.

The kiln does not reach the programmed temperatures. The controller displays an error message. Check the heating elements for visible damage, e.g. cracks. Heating elements might be broken and therefore not working properly!

Have a qualified electrician check if all heating elements are properly connected.

8.2. The controller displays an event message

You will find the relevant explanation in your user's manual for the controller.

9. WARRANTY OBLIGATIONS

In accordance with the legal warranty regulations we provide a 24-month warranty from date of invoice.

As well as the heating elements (subject to wear) the following are excluded from the scope of warranty:

- Damage caused by the customer.
- Damage caused by the fired material, e.g. due to temperature limits being exceeded.
- • Corrosion caused by aggressive glazes or insufficient ventilation of the firing chamber.
- Damage caused by improper transport.
- Damage due to chemical reactions during firing for which the kiln is not intended (such as salt glaze).
- The manufacturer is not liable for any damage resulting from improper operation.

Please note: The firebricks of the kiln lining are exposed to significant temperature fluctuations. This may cause hairline cracks in the firebrick lining. This is normal and does not affect the functionality of the kiln. It cannot therefore be accepted as a reason for complaint.

What to do in the case of warranty/damage:

Please notify us - before any costs arise. If any claims arise, please state the kiln type, product number and the date of purchase or the year of construction (see type plate on switch cabinet).

PLEASE NOTE: We refer to the General Terms and Conditions (as of 04/2014) of Helmut Rohde GmbH.

10. PROPERTY RIGHTS / TRADE NAMES / DISCLAIMER

The contents of the instruction manual are purely informative. Changes may be made without prior notice.

We mention names, trade names, product identifications etc. without special identification, as they are generally known. Those names and identifications, however, may be the property of companies or institutions and subject to copyright.

11. CONTACTS / ASSISTANCE

If you have any questions regarding your frontloader melting kiln, spare parts and additional equipment, please contact your local dealer.

Enjoy working with your new kiln! We wish you excellent firing results.
Your kiln team